

Qualys Gateway Service v3.x

Release Notes

Version 3.1.0

November 30, 2022

Qualys Gateway Service 3.1 brings you the following enhancements.

Visibility to Appliance Upgrade Details and Primary Disk and RAM details

With this release, you can now view whether the appliance is upgraded to the latest services or not. From the appliance's Text User Interface (TUI), you can click the **Info tab** and click **OK** to see the details.



REPOSITORY	TAG	IMAGE ID
cams-logstash	1.5.3-14	c50aa1f08
cams-rsyslog	1.5.3-14	52ec18f28
cams-metrics	1.5.0-43	6f74eff22
cams-haproxy	1.5.0-43	a170ff2f8
cams-squid	1.5.0-43	c82b3908e

↓ (+)

< EXIT >

Also, when the minimum requirement for the primary disk and RAM are not fulfilled, the following message is shown on the appliance TUI under the Info tab.

Note: The minimum requirement for the primary disk is 30 GB, and the RAM is 16 GB. You can click the **Info tab** and click **OK** to see the message.



```
Requirements mismatch :
This appliance does not meet the minimal requirements and
may not function correctly.

WARNING : Total disk space is 20 GB. Minimum disk space
requirement of 30 GB is not matched
WARNING : Found 8.1709 GB installed memory. Minimum RAM
requirement of 16 GB is not matched

Please refer to the CAMS setup instructions

< OK >
```

Additionally, you can now see the following image versions on the appliance TUI under the Info tab.

Go to **Diagnostics > Images and Diagnostics** and click **OK** to see the image versions.

REPOSITORY	TAG	IMAGE ID	CREATED	SIZE
cams-logstash	1.5.3-14	c50aa1f08069	4 weeks ago	759MB
cams-rsyslog	1.5.3-14	52ec18f28fca	4 weeks ago	84.4MB
camsd	1.5.3-14	a04f16fa7e24	4 weeks ago	794MB
cams-metrics	1.5.0-43	6f74eff22f12	5 months ago	95.7MB
cams-haproxy	1.5.0-43	a170ff2f846c	5 months ago	87.6MB

✔✔✔

74%

< EXIT >